



A new contract research service expands the in vitro ADME assays that can be analyzed using RapidFire technology. The extension of the existing CYP Inhibition service now includes analysis of a range of assays including microsomal stability, permeability, and plasma protein binding. Additionally, custom services are offered for the development of analysis methods for other ADME assays. Clients can preserve their assay method integrity as assays are processed using their own laboratory methods and then analyzed using RapidFire mass spectrometry technology. RapidFire processes samples at speeds of up to six seconds per sample, enabling high throughput and good data turnaround.

Company: Biocius Life Sciences