



The KingFisher Duo nucleic acid purification system is a low- to medium-throughput package, which includes the KingFisher magnetic particle processor instrument, optimized DNA/RNA isolation kits, plastics consumables, and BindIt software to provide a fully automated workflow. The KingFisher Duo has a small benchtop footprint, making it useful for space-restricted laboratories. With the ability to process a variety of sample types in volumes up to 5 mL, users can obtain high yields of pure product for downstream analysis. The system is suitable for research and diagnostic labs that require nucleic acid purification, and it can be modified for use in proteomics and cell isolation.

Company: Thermo Fisher Scientific